**I ♥ ZEBRA
SUPPLIES**

LEOSS
moč podatkov

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Suggested Ribbons for Applications with Chemicals

	Weak					Moderate				Harsh			Extreme				
	Blood	Body Fluid	Salt Water	Water	Window Cleaner	Alcohol	Ammonia	Bleach	IPA	Gasoline	Grease	Oil	Acetone	IR Reflow	MEK	TCE	Xylene
5095	X	X	X	X	X		X	x	X								

"X" indicates acceptable chemical resistance

**I ♥ ZEBRA
SUPPLIES**

LEOSS
moč podatkov

Cryogenic Testing: Test Procedure

Labels were applied to glass vials (2.8 cm OD), polypropylene centrifuge tubes (3.5 cm OD, 50ml) and glass microscope slides and allowed a 24 hour dwell time before exposure to above conditions.

Environment	Test Method	Typical Results
High Temp.	30 days at listed temperature	No visible effect at 90 C (194 F)
Low Temp.	30 days at -70 C (-94 F)	No Visible effect
Freezer	3 cycles of 16 hours at -70 C (-94 F)/ 8 hours at room temp.	Glass vial : Recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Pressure Cooker	3 cycles of 1 hour in 121 C (250 F) 15 psi pressure cooker/ 23 hours room temperature	Glass vial : Recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Liquid Nitrogen	3 cycles of 4 hours at - 196 C (-320 F)/ 20 hours at room temperature	Glass vial : Not recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Freezer to boiling water	1 hour at -70 C (-94 F) then placed in boiling water 100 C (212 F)	Glass vial : May work, must test PP centrifuge tube: Recommended Glass microscope slide: May work, must test Flat PP: Recommended
Liquid Nitrogen to boiling water	1 hour at -196 C (-320 F) then placed in boiling water 100C (212 F) for 10 minutes	Glass vial : Not Recommended PP centrifuge tube: Recommended Glass microscope slide: May work, must test Flat PP: Recommended

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